

MEDSOURCE DESKTOP STEAM STERILIZER

YLT-250A

YLT-255A



ISO9001, ISO13485, and GMP.

www.medsource.com.tw



MEDSOURCE CO., LTD.

MEDSOURCE DESKTOP STEAM STERILIZER

Function Features:

1. The housing of the sterilizer is composed of upper lid coated with baking enamel, and lower lid covered by stainless steel for protection. It has advantages of corrosion resistance, cleaning convenience and aesthetics, etc.
2. It's an automatic controlled process including heating up, sterilization, steam exhaust and drying, and its operation process will be displayed by the control panel.
3. Control parts including: temperature adjuster, sterilization timer, drying timer, pressure gauge, water inlet and water drainage device, these can all be operated and monitored at the front of the stand.
4. The integrated door cover is made of PC synthetic plastic materials, with advanced and aesthetic spiral T-handler made of special steel. It has the advantages of simple and convenient operation, high safety, meet ergonomic principles, and long life cycle etc.
5. Exhaust valve of heavy mallet type can make sterilization and dry function more perfect.
6. The international standard design of safety valve will take action immediately when the pressure of inner chamber exceeds the 10% of max. usage pressure value, which meets the pressure capacity safety regulation.
7. Quality standard: certification of ISO9001, ISO13485, and GMP.

Specification

Functions	Model	YLT-250A	YLT-255A
Exterior Dimension (cm)		(W)47x(D)61X(H)39cm	(W)47x(D)73X(H)39cm
Inner Chamber Dimension (cm)		(Ø)25x (D)41cm	(Ø)25x (D)50cm
Inner Chamber Capacity (L)		20L	24L
Consumption Power (W)		Inner Chamber Heater 1450W, Dry Heater Piece 300W	Inner Chamber Heater 1800W, Dry Heater Piece 300W
Supply Voltage		230VAC 1phase 8Amp	230VAC 1phase 10Amp
Water Tank Capacity (L)		5L	5L
Net Weight (KG)		39KG	42KG
Working Temperature		121~134°C	121~134°C
Operate Pressure		Max 2.1Kg/cm ²	Max 2.1Kg/cm ²
Standard Equipment		1 of Stainless Steel Tray	1 of Stainless Steel Tray
Optional Accessories		<ul style="list-style-type: none"> • Tray Frame x one pc • Tray x 4 pcs (2 pcs for 375x210x20mm, 2 pcs for 375x165x20mm) • Tray holder x one pc 	<ul style="list-style-type: none"> • Tray Frame x one pc • Tray x 4 pcs (2 pcs for 475x210x20mm, 2 pcs for 475x165x20mm) • Tray holder x one pc



MEDSOURCE CO., LTD.

4F., No.3, Ningbo W. St., Zhongzheng Dist., Taipei City 100, Taiwan (R.O.C.)

TEL: +886-2-23930557

FAX: +886-2-23930557 • 886-2-23932130

E-mail: medsource@medsource.com.tw medsource1688@medsource.com.tw

www.medsource.com.tw