


E.N.T Treatment Table & Chair



Electric treatment chair for ENT

MODEL NO.	YL-6500
OVER ALL	170(L) x 63(W) x 58-78(H) cm
BACK SECTION	90° fixed
MATTRESS	PVC foam and cover with PVC plastic
ELEVATION	Approx. 58cm to 78cm adjustable by elec. Hydraulic pump
HEADREST	Adjusted by manual control
POWER	110V or 220V or 240V 50/60Hz
CONSUMPTION	110V / 7A 220V / 4A 240V / 4A
BASIC BASE DIA.	23"

E.N.T Treatment Table Unit

MODEL NO.	YL-7600S
OVER ALL	70(L) X 45(W) X 85(H) CM
HOUSING	COVER WITH ST.ST.#304
SPRAY GUN	3 PCS.
AIR COMPRESSOR	OIL-LESS PUMP
SUCTION PUMP	OIL-LESS PUMP
LIGHT	MOVING FREE ARM LIGHT (100W)
MEDICAL CUP	8 PCS
COLLECTORS (ST.#304)	2PCS
POWER CONSUMPTION	110V/10A, 220V/5A, 240V/5A
POWER SOURCE	110V,OR 220V, OR240V 50/60HZ



MEDSOURCE CO., LTD.

4F., No.3, Ningbo W. St., Zhongzheng Dist., Taipei City 100, Taiwan (R.O.C.)

TEL: +886-2-23930557 • 886-2-23935000

FAX: +886-2-23930557 • 886-2-23932130

E-mail: medsource@medsource.com.tw

medsource1688@medsource.com.tw

www.medsource.com.tw